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Revision changes

26-03-2019 biobank tracking system is changed in Biobank management system and a comment has been added.

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1. Goal

To describe how serum 10ml must be prepared and stored for purposes of the Radboud Biobank, so that:

1. Interested parties including potential users know how the biomaterial has been handled.
2. The laboratory can assess whether they can process the biomaterial as described under paragraph 2 "Preparation".
3. The laboratory knows the requirements for the collection, processing and storage of the biomaterial to enable registration of deviations.
4. The sub biobank knows the requirements for the collection and transport of the biomaterial to ensure correct delivery.
5. The RB is able to couple this procedure to the biomaterial in storage to determine the fitness for purpose.

1.1 Scope of application

This procedure is applicable to all employees concerned at the Radboud Biobank, the Sub biobank/Project and the the Radboud Laboratory for Diagnostics.

2. Protocol serum 10 ml

Collection

Delivery	In collection tube.
Type of tube	Standard tube 10 ml, without gel, without protease inhibitors. BD-tube 367896 or equivalent.
Temperature	Keep tubes at room temperature until preparation for storage.

Preparation

Time until freezing	Prepare and freeze the collected blood as soon as possible. Goal: within 2 hours. Maximum: within 4 hours.
Coagulation time	Minimum coagulation time is 60 minutes, maximum 120 minutes. Coagulation takes place at room temperature.
Centrifuge	g-force: 2000g Time: 10 minutes. Temperature: At room temperature.

Storage

Register	Register the Microtubes in the biobank management system and attach a cryovial label.
Aliquotation	Store the biomaterial in 0.5 ml aliquots to a maximum of 6 aliquots. If there is biomaterial left over add it to the 6 aliquots.

Type of tube 2 ml polypropylene microtube and screw-top with O-ring. GREI722301UMC and GREI368380UMC or equivalent.

Temperature -80°C

Miscellaneous

Comments The MEN pearl uses a SST tube with silica clot activator, 8.5 ml, BD-tube 367988 or equivalent.

Deviations Note all deviations from this procedure in the biobank management system. Assess the samples macroscopically for the following deviations and register them in the Biobank management system.

- Haemolytic
- Icteric
- Lipemic

3. Quality assurance

Review takes place via i) an internal audit (see the [audit year plan](#) in DMS) and ii) 2-monthly performance reviews of deviations registered in the biobank management system (DMS map "[Registraties en controles](#)").

3.1 Performance indicator

95% of the sub biobank samples have been prepared and stored within the parameters of this procedure.

4. Comments

-

5. Relevant documents

Qdocs, literature, legislation and codes of conduct, website etc.

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