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## Revision changes

### 1. Goal

To describe how DNA from umbilical cord blood must be prepared and stored for purposes of the Radboud Biobank (RB), so that:

1. Interested parties including potential users know how the biomaterial has been handled.
2. The laboratory can assess whether they can process the biomaterial as described under paragraph 2 "Preparation".
3. The laboratory knows the requirements for the collection, processing and storage of the biomaterial to enable registration of deviations.
4. The sub biobank or project knows the requirements for the collection and transport of the biomaterial to ensure correct delivery.
5. The RB is able to couple this procedure to the biomaterial in storage to determine the fitness for purpose.

#### 1.1 Scope of application

This procedure is applicable to all employees concerned at the Radboud Biobank, the Sub biobank/Project and the Department of Human Genetics of the Radboudumc.

## 2. Protocol DNA isolation from umbilical cord blood

### Collection

**Collection** Collect blood from umbilical cord (ex utero) as fast as possible after the placenta has detached before the clotting process begins. The process of taking umbilical cord blood is described in Qdoc 015297. For DNA isolation it does not matter whether the blood comes from vein, artery or combination of both.

**Delivery** In collection tube.

**Type of tube** Standard EDTA tube, no gel and protease inhibitors.  
At least 4 ml; 6ml is preferred.

**Type of needle** Butterfly needle or use UmbiliCup (DeRoyal, 72-8000).  
UmbiliCup is a (safety) umbilical cord blood collection system used for umbilical cord blood sampling without utilizing an exposed sharp needle.

**Temperature** Keep at room temperature (RT) until isolation or store at -80°C when isolation is not possible within 48 hours.

### Preparation

**Time until freezing** Isolate DNA within 48 hours of collection, if blood is stored at -80 °C isolate DNA within a maximum of 4 months.

**Isolation protocol** DNA isolation from EDTA-whole blood with ChemagicStar (Qdoc 047239).

**DNA-concentration and volume** After isolation determine the DNA-concentration (in ng/μl) and the OD-ratio 260-280 nm from the stock solution and register in the Biobank management system.

### Storage

**Aliquotation** Store the DNA in a 500 μl stock solution (1 aliquot).  
N.B. The first time the DNA is issued make a work solution of 200 μl and 100 ng/μl and store in MATRIX-tube.

**Title: SOP DNA from umbilical cord blood**

Type of tube                    MATRIX-tube (Matrix™ 2D Barcoded Open-Top Storage Tubes, ThermoFisher Scientific).

Temperature                    -20°C

**Miscellaneous**

Data management            Register the following data in the biobank management system:  
 1. Date and time blood was collected  
 2. Input-date is time of storage at -20°C after DNA-isolation.  
 3. Identification number (pseudomised)  
 4. Name of sub-biobank/ pearl/ project.  
 5. OD-ratio and DNA-concentration.  
 6. Date and time of work solution is stored at -20°C (is the date the solution is made), volume and concentration of the working solution.  
 7. Location of the DNA in the storage facility.  
 8. Origin DNA sample (Blood or Saliva).

Deviations                    Note all deviations from this procedure in the Biobank management system

**3. Quality assurances**

Review takes place via i) an internal audit (see the [audit year plan](#) in DMS) and ii) 2-monthly performance reviews of deviations registered in the biobank management system (DMS map "[Registraties en controles](#)").

**3.1 Performance indicator**

95% of the sub biobank samples have been prepared and stored within the parameters of this procedure.

**4. Comments**

-

**5. Relevant documents**

Qdocs, literature, legislation and codes of conduct, website etc.

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