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Revision changes

1. Goal

To describe how FFPE and frozen tissue must be prepared and stored within in the Radboud Biobank, so that:

1. Potential users know how the biomaterial has been handled.
2. The laboratory preparing the biomaterial for storage knows the minimal requirements for the handling of it.

1.1 Scope of application

This procedure is applicable to all employees concerned at the Radboud Biobank and the Radboud Laboratory for Diagnostics and the Pathology department.

2. Protocol FFPE and frozen tissue

Collection

| | |
|-------------------|---|
| Delivery | Fresh tissue without fixative |
| Type of container | No specific requirements only size of the container must fit with the size of the tissue. |
| Temperature | Store tissue at room temperature until preparation. |

Preparation

| | |
|-----------------------|--|
| Time until processing | As quickly as possible. |
| Temperature | Room Temperature. |
| Preparation protocol | Tissue is assessed, if there is sufficient material: <u>Diagnostic residual material:</u> Minimal collection: 2x (not affected) normal tissue (1x frozen / 1x paraffin) 2x (tumor) not normal tissue (1x frozen / 1x paraffin) Preferred collection as much as possible (max 16): 8x (not affected) normal tissue (4x frozen / 4x paraffin) 8x (tumor) not normal tissue (4x frozen / 4x paraffin) <u>Non-diagnostic residual material:</u> It is sub Biobank dependent what and how much is taken. All the material is processed in separate tubes. |

Storage

| | |
|--------------|--|
| Register | Register the material in the Biobank management system and attach a cryovial label on the frozen tissue tubes; for paraffin tissue transfer the item ID on the holder. |
| Aliquotation | Not applicable. |
| Type of tube | Frozen tissue: 1,8 ml Cryotube (Thermo Scientific) VWR 479-6837 or equivalent. Paraffine tissue: Not applicable. |

Temperature Frozen tissue: -80°C
 Paraffin tissue: Room temperature

Miscellaneous

Deviations Note all deviations from this procedure in the Biobank management system.

3. Quality assurance

Review takes place via an internal audit (see the audit plan).

3.1 Performance indicator

95% of the relevant samples of the sub Biobank have been prepared and stored within the parameters of this procedure.

4. Comments

5. Relevant documents

Qdocs, literature, legislation and codes of conduct, website etc.

| link | Titel |
|---------------------------|---|
| Qdoc 065989 | Processchema team 2 |
| Qdoc 075999 | Processchema Biobank |
| Qdoc 065547 | Logboek Biobank aanvragen |
| Bio-16206 | Monstermanagement RLD t.b.v. Radboud Biobank |
| Bio-17159 | Monstermanagement Pathologie t.b.v. Radboud Biobank |

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