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Revision changes

1. Goal

To describe how 24h urine must be prepared and stored for purposes of the Radboud Biobank, so that:

1. Interested parties including potential users know how the biomaterial has been handled.
2. The laboratory can assess whether they can process the biomaterial as described under paragraph 2 "Preparation".
3. The laboratory knows the requirements for the collection, processing and storage of the biomaterial to enable registration of deviations.
4. The sub biobank knows the requirements for the collection and transport of the biomaterial to ensure correct delivery.
5. The RB is able to couple this procedure to the biomaterial in storage to determine the fitness for purpose.

1.1 Scope of application

This procedure is applicable to all employees concerned at the Radboud Biobank, the Sub-biobank/Project and the Radboud Laboratory for diagnostics.

2. Protocol 24h urine

Collection

Collection

- Discard the first morning urine.
- Note the date and time of start collection.
- Collect all urine over the next 24 hours including the last specimen at the same time (start time) the next day.
- Urine will be saved and collected in the container(s) provided

Transport

Transport the urine to the Radboudumc as soon as possible after finishing the collection.

Aim: The same day the collection is finished.

Type of tube

3 liter urine container, VWR 44333-016 or equivalent.

Temperature

Keep urine during collection at 2-8 °C in the fridge.

Preparation

Time until freezing

Prepare and freeze the urine as soon as possible.

Aim: within 2 hours

Maximum: within 4 hours after receiving the urine on the Radboudumc.

Sample handling

Determine the total weight of the urine and register in GLIMS.

Mix the urine and fill out in centrifugation tubes, Greiner 455001 or equivalent.

Centrifuge

g-force: 2000 g.

Time: 10 minutes.

Temperature: At Room Temperature

Storage

Register

Register the 2 ml (micro) tubes in the Biobank management system and attach a cryovial label.

Title: SOP 24h Urine

Radboud Biobank

Aliquotation	Store the urine in 6 aliquots of minimal 0,9 ml (and discard the remaining material).
Type of tube	2 ml polypropylene microtube with screw-top with O-ring. GREI722301UMC and GREI368380UMC or equivalent.
Temperature	-80 °C.

Miscellaneous

Deviations	Note all deviations from this procedure in the Biobank management system. Give a summary of the deviations important for this description. <ul style="list-style-type: none"> • Urine collection is not complete. • Urine is not cooled during collection
------------	---

3. Quality assurances

Review takes place via i) an internal audit and ii) 2-monthly performance reviews of deviations registered in the biobank management system.

3.1 Performance indicator

95% of the sub biobank samples have been prepared and stored within the parameters of this procedure.

4. Comments

-

5. Relevant documents

Qdocs, literature, legislation and codes of conduct, website etc.

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