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07-07-2019 Title is changed

07-07-2019 Performance indicator is changed

26-03-2019 Time till preparation is changed in time until freezing.

26-03-2019 biobank tracking system is changed in Biobank management system

## 1. Goal

To describe how urine must be prepared and stored within the framework of the Parelinoer Institute for the purposes of the Radboud Biobank, so that:

1. Potential users know how the biomaterial has been handled.
2. The laboratory preparing the biomaterial for storage knows the minimal requirements for the handling of it.

### 1.1 Scope of application

This procedure is applicable to all employees concerned at the Radboud Biobank and the Radboud Laboratory for Diagnostics.

## 2. Protocol urine standard

### Collection

Delivery In urine jar or large polypropylene tube.

Temperature Keep on ice or at 4°C until preparation.

### Preparation

Time until freezing Maximum of 4 hours after collection.

Centrifuge Mix sample before processing.  
g-force: 2000g  
Time: 10 minutes.  
Temperature: 4°C.

### Storage

Register Register the microtubes in the biobank management system and attach a cryovial label.

Aliquotation Store the biomaterial in 0.9 ml aliquots to a maximum of 6 aliquots. If there is biomaterial left over add it to the 6 aliquots.

Type of tube 2 ml polypropylene microtube and screw-top with O-ring.  
GREI722301UMC and GREI368380UMC or equivalent.

Temperature -80°C

### Miscellaneous

Deviations Note all deviations from this procedure in the biobank management system.

## 3. Quality assurance

Review takes place via an internal audit (see the audit plan).

### 3.1 Performance indicator

95% of sub Biobank samples have been prepared and stored within the parameters of this procedure.

## 4. Comments

-

## 5. Relevant documents

Qdocs, literature, legislation and codes of conduct, website etc.

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