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1. Goal

To describe how feces must be prepared and stored for purposes of the Radboud Biobank (RB), so that:

1. Interested parties including potential users know how the biomaterial has been handled.
2. The laboratory can assess whether they can process the biomaterial as described under paragraph 2 "Preparation".
3. The laboratory knows the requirements for the collection, processing and storage of the biomaterial to enable registration of deviations.
4. The sub biobank or project knows the requirements for the collection and transport of the biomaterial to ensure correct delivery.
5. The RB is able to couple this procedure to the biomaterial in storage to determine the fitness for purpose.

1.1 Scope of application

This procedure is applicable to all employees concerned at the Radboud Biobank, the Sub biobank/Project and the Radboud laboratory for diagnostics.

2. Protocol Feces collection in tube with spoon

Collection

Delivery	In feces tube with spoon.
Type of tube	Sterile specimen container with spoon, 20 X 76 mm with polypropylene screw cap (Sarstedt 80.734.001).
Instructions	Spoon 1 gram of feces (approx. 1ml) - clean collected – in the feces collection tube.
Temperature	When feces is collected at home keep the tube in the refrigerator and hand it in at the Radboudumc as soon as possible.

Preparation

Time until freezing	Freeze the collected feces as soon as possible. Aim: within 4 hours (collection in Radboudumc) Maximum: within 72 hours
---------------------	---

Storage

Register	Register the feces tubes in the biobank management system and attach a cryovial label.
Type of tube	Sterile specimen container with spoon, 20 X 76 mm with polypropylene screw cap (Sarstedt 80.734.001).
Temperature	-80 °C.

Miscellaneous

Deviations Note all deviations from this procedure in the Biobank management system.

3. Quality assurances

Review takes place via i) an internal audit (see the [audit year plan](#) in DMS) and ii) 2-monthly performance reviews of deviations registered in the biobank management system (DMS map "[Registraties en controles](#)").

3.1 Performance indicator

95% of the sub biobank samples have been prepared and stored within the parameters of this procedure.

4. Comments

-

5. Relevant documents

Qdocs, literature, legislation and codes of conduct, website etc.

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