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Revision changes

Translation oversight corrected

26-03-2019 biobank tracking system is changed in biobank management system

1. Goal

To describe how EDTA whole blood for DNA extraction must be prepared and stored for the purposes of the Radboud Biobank, so that:

1. Potential users know how the biomaterial has been handled.
2. The laboratory preparing the biomaterial for storage knows the minimal requirements for the handling of it.

1.1 Scope of application

This procedure is applicable to all employees concerned at the Radboud Biobank and the Radboud Laboratory for Diagnostics.

2. Protocol EDTA whole blood for DNA extraction

Collection

Delivery In collection tube.

Type of tube Standard 6 ml tube EDTA K2E, without gel, without protease inhibitors.
BD tube 367864 or equivalent.

Temperature Keep tubes at room temperature until prepared for storage.

Preparation

Time until freezing Maximum: 96 hours.

Centrifuge Not applicable.

Storage

Registration Register the EDTA whole blood tubes in the biobank management system and attach a cryovial label.

Aliquotation Store the complete EDTA whole blood tube without further preparation.

Temperature -80°C

Miscellaneous

Deviations Note all deviations from this procedure in the biobank management system.

3. Quality assurance

Review takes place via an internal audit (see the audit plan).

3.1 Performance indicator

95% of the relevant samples of the sub biobank have been prepared and stored within the parameters of this procedure.

4. Comments

-

5. Relevant documents

Qdocs, literature, legislation and codes of conduct, website etc.

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