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## Revision changes

Concentration citrate is added.

25-03-2018: Biobank tracking system is changed to Biobank management system

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## 1. Goal

To describe how citrate plasma (9 ml) must be prepared and stored for purposes of the Radboud Biobank, so that:

1. Interested parties including potential users know how the biomaterial has been handled.
2. The laboratory can assess whether they can process the biomaterial as described under paragraph 2 "Preparation".
3. The laboratory knows the requirements for the collection, processing and storage of the biomaterial to enable registration of deviations.
4. The sub biobank knows the requirements for the collection and transport of the biomaterial to ensure correct delivery.
5. The RB is able to couple this procedure to the biomaterial in storage to determine the fitness for purpose.

### 1.1 Scope of application

This procedure is applicable to all employees concerned at the Radboud Biobank, the Sub biobank/Project and the Radboud Laboratory for Diagnostics.

## 2. Protocol citrate plasma 9 ml

### Collection

Fasting patient	Dependant on the requirements of the sub-biobank/Project.
Delivery	In collection tube.
Type of tube	Standard 9 ml tube, 3,2% sodium citrate, no glass tubes unless silicone lined, without gel, without protease inhibitors. Greiner tube 455322 or equivalent.
Temperature	Keep tubes at room temperature until prepared for storage.

### Preparation

Time until freezing	Prepare and freeze the collected blood as soon as possible.. Aim: within 2 hours. Maximum: within 4 hours.
Centrifuge	g-force: 2000g Time: 10 minutes. Temperature: At room temperature.

### Storage

Register	Register the microtubes in the Biobank management system and
----------	--

attach a cryovial label.

Aliquotation            Store the biomaterial in 0.5 ml aliquots to a maximum of 6 aliquots. If there is biomaterial left over add it to the 6 aliquots.

Type of tube            2 ml polypropylene microtube and screw-top with O-ring. GREI722301UMC and GREI368380UMC or equivalent.

Temperature            -80°C

**Miscellaneous**

Deviations            Note all deviations from this procedure in the Biobank management system.

**3. Quality assurance**

Review takes place via i) an internal audit (see the [audit year plan](#) in DMS) and ii) 2-monthly performance reviews of deviations registered in the biobank management system (DMS map [“Registraties en controles”](#)).

**3.1 Performance indicator**

95% of the sub biobank samples have been prepared and stored within the parameters of this procedure.

**4. Comments**

-

**5. Relevant documents**

Qdocs, literature, legislation and codes of conduct, website etc.

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