Title: SOP EDTA-plasma in cryovial

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Revision changes

07-01-2020 Performance indicator is changed
 07-01-2020 Registration of deviations is changed
 30-11-2021 Format is changed
 30-11-2021 Relevant documents are updated

1. Goal

To describe how EDTA plasma in cryovial must be prepared and stored for purposes of the Radboud Biobank, so that:

1. Interested parties including potential users know how the biomaterial has been handled.
2. The laboratory can assess whether they can process the biomaterial as described under paragraph 2 “Preparation”.
3. The laboratory knows the requirements for the collection, processing and storage of the biomaterial to enable registration of deviations.
4. The sub biobank knows the requirements for the collection and transport of the biomaterial to ensure correct delivery.
5. The RB is able to couple this procedure to the biomaterial in storage to determine the fitness for purpose.

1.1 Scope of application

This procedure is applicable to all employees concerned at the Radboud Biobank and the sub-biobank/project.

2. Protocol EDTA-plasma in cryovial

Collection

Delivery In 2ml microtube (Batch wise delivery).

Preparation

Time until freezing Prepare and freeze the plasma on site as soon as possible.
 Aim: within 2.hours after collection.
 Maximum: within 4 hours.

Temperature -80°C

Transport On dry ice.

Storage

Registration Register the microtubes in the biobank management system.

Type of tube 2 ml polypropylene microtube and screw-top with O-ring.
 GREI722301UMC and GREI368380UMC or equivalent.

Temperature -80°C

Miscellaneous

Deviations Note all deviations from this procedure in the biobank management system.
 (give a summary of the deviations important for this description)

- Haemolytic
- Icteric
- Lipemic

Title: SOP EDTA-plasma in cryovial**3. Quality assurance**

Review takes place via i) an internal audit and ii) 2-monthly performance reviews of deviations registered in the biobank management system.

3.1 Performance indicator

95% of the relevant samples from the sub-biobank have been prepared and stored within the parameters of this procedure.

4. Comments

-

5. Relevant documents

Qdocs, literature, legislation and codes of conduct, website etc.

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