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## Revision changes

### 1. Goal

To describe how a SARS-CoV-2 throat/nasopharynx swab must be collected and stored for Radboud Biobank purposes so that:

1. Interested parties including potential users know how the biomaterial has been handled.
2. The laboratory can assess whether they can process the biomaterial as described under paragraph 2 "Preparation".
3. The laboratory knows the requirements for the collection, processing and storage of the biomaterial to enable registration of deviations.
4. The sub biobank knows the requirements for the collection and transport of the biomaterial to ensure correct delivery.
5. The RB is able to couple this procedure to the biomaterial in storage to determine the fitness for purpose.

#### 1.1 Scope of application

This procedure is applicable to all employees concerned at the Radboud Biobank, the sub-biobank/project and the Radboud Laboratory for Diagnostics.

## 2. Protocol swab collection

### Collection

Collection	See collection instruction in <a href="#">Qdoc 080291</a> . One cotton swab is used for a deep throat swab. Place the cotton swab in the supplied transport tube and cut or break it off. A second swab is used for a deep nasopharynx swab. This swab is placed in the <b>same</b> transport tube as the throat swab, cut or break it off and close the lid.
Delivery	In collection tube.
Type of tube	3 ml UTM (Universal Transport Medium) nasopharyngeal and regular FLOQ swab, 321 C
Temperature	Room temperature.

### Preparation

Time until freezing	Preferably freeze the collected swab as soon as possible. Maximum: within 4 hours after collection.
---------------------	--

### Storage

Register	Register the collection tube in the Biobank management system and attach a cryovial label. Remove first all labels from the swab tube before adding the cryovial label
Aliquotation	Do not divide into aliquots but freeze the complete tube in an upright position.
Temperature	-80°C.

### Miscellaneous

Deviations	Note all deviations from this procedure in the Biobank management system.
------------	---

### 3. Quality assurance

Review takes place via i) an internal audit and ii) 2-monthly performance reviews of deviations registered in the biobank management system.

#### 3.1 Performance indicator

95% of the sub biobank samples have been prepared and stored within the parameters of this procedure.

### 4. Comments

-

### 5. Relevant documents

Qdocs, literature, legislation and codes of conduct, website etc.

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